

ad notam POWER ADAPTOR

ET_523

External power adaptor for **ad notam** monitor systems.

Applicable with **ad notam** article numbers:

FPD-0216




PRODUCT FEATURES

- Universal AC Input
- Meets CEC & ErP Level V
- Meets EN55022 Class "B"
- Meets CISPR/FCC Class "B"
- Continuous Short Circuit Protection
- High Efficiency Operation, 87% Typical
- Over Voltage Protection

PACKAGE INCLUDES

- 1 x power adaptor
- 1 x power cord

Power	Voltage Input	Voltage Output	Operating Temperature	Humidity	Weight
45W / 1.88A	100 ~ 240V AC	24V DC	0 - 40 °C	20 - 90%	250g



ARTICLE NO.

Art.No.	Code
ET_	523